
# Curve fitting in HyperSpy: Applications to EELS data analysis

**Josh Taillon and Andy Herzing** 

April 5, 2018





### Create some test data

• 50x50 array of Gaussians with randomly varying A, mu, and sigma


```
In [2]: amps = hs.signals.Signal2D(np.random.randint(500,1500,size=[20,20]))
    centers = hs.signals.Signal2D(np.random.randint(45,55,size=[20,20]))
    sigmas = hs.signals.Signal2D(np.random.randint(5,10,size=[20,20]))


    xaxis = np.tile(np.arange(0,100),20**2)
    xaxis = xaxis.reshape([20,20,100])
    xaxis = np.rollaxis(xaxis,2)
    curves = amps.data/np.sqrt(2*np.pi*sigmas.data**2)*np.exp(-(xaxis-centers.data)**2/(2*sigmas.data**2))
    curves = np.rollaxis(curves,0,3)
    gaussians = hs.signals.SignallD(curves)

gaussians
```



In [3]: gaussians.plot()








### Create a model


### Add single Gaussian component:



### Perform fit at all pixel locations:


```
In [6]: gaussians_model.multifit()
    gaussians_model.plot()
```









## Compare fitted values to test data









### Add some noise to the test data

```
In [8]: noisy = gaussians.deepcopy()
noisy.unfold()

for i in range(0,20**2):
    noisy.data[i,:] =noisy.data[i,:] + np.random.poisson(np.sqrt(noisy.data[i,:].max()),100)
noisy.fold()
```



# Modeling the noisy data

Create model and add a single Gaussian component:

```
In [9]: noisy_model = noisy.create_model()
noisy_model.append(hs.model.components1D.Gaussian())
```

Add a linear offset to account for 'background' from adding noise:

```
In [10]: noisy_model.append(hs.model.components1D.Offset())
```

Estimate model parameters using a quick arithmetic estimation:


```
In [11]: noisy_model['Gaussian'].estimate_parameters(signal=noisy, x1=10, x2=90, only_current=False)
noisy_model['Offset'].estimate_parameters(signal=noisy, x1=10, x2=90, only_current=False);
```




### Perform fit at all image pixels:

```
In [12]: noisy_model.multifit(show_progressbar=True)
    noisy_model.plot(plot_components=True)
```

100% 400/400 [00:03<00:00, 124.82it/s]







# More Realistic Application:

Fitting EELS Fine Structure



### **Useful aside: The EELS Database**



- <a href="https://eelsdb.eu/">https://eelsdb.eu/</a>
  - Largest open-access repository of EELS and X-ray absorption spectra
  - Nearly 300 spectra covering 43 elements and growing
- Database can be directly queried by:
  - title "Hexagonal Boron Nitride"
  - formula "BN"
  - $\blacksquare$  element "B,N"
  - edge "K" or "L1"
  - Etc., many other options



### Spectra can be pulled directly from database into HyperSpy (requires internet connection):



```
In [14]: # Search my element and spectrum type:
         hs.datasets.eelsdb(element=('Fe','O'), spectrum type='coreloss')
Out[14]: [<EELSSpectrum, title: Iron (III) oxide, hematite, dimensions: (|1024)>,
          <EELSSpectrum, title: Iron Oxide, dimensions: (|1024)>,
          <EELSSpectrum, title: Hematite, dimensions: (|1024)>,
          <EELSSpectrum, title: Iron Oxide Hematite, dimensions: (|926)>,
          <EELSSpectrum, title: Iron Oxide Hematite, dimensions: (|620)>,
          <EELSSpectrum, title: Ti-ferrite (spinel), dimensions: (|780)>,
          <EELSSpectrum, title: Ti-ferrite (spinel), dimensions: (|780)>,
          <EELSSpectrum, title: Iron Oxide Magnetite, dimensions: (|911)>,
          <EELSSpectrum, title: Iron Oxide Magnetite, dimensions: (|629)>,
          <EELSSpectrum, title: Iron Oxide Siderite, dimensions: (|826)>,
          <EELSSpectrum, title: Iron Oxide Siderite, dimensions: (|605)>,
          <EELSSpectrum, title: Iron Oxide 2-lines ferrihydrite, dimensions: (|671)>,
          <EELSSpectrum, title: Iron Oxide 2-lines ferrihydrite, dimensions: (|660)>,
          <EELSSpectrum, title: Iron Oxide Goethite, dimensions: (|1024)>,
          <EELSSpectrum, title: Iron Oxide Goethite, dimensions: (|704)>,
          <EELSSpectrum, title: Iron Oxide Cl-containing akaganeite, dimensions: (|427)>,
          <EELSSpectrum, title: Iron Oxide Cl-containing akaganeite, dimensions: (|552)>,
          <EELSSpectrum, title: ferrous Titanate, dimensions: (|1024)>,
          <EELSSpectrum, title: Iron titanium oxide, dimensions: (|1024)>,
          <EELSSpectrum, title: LaFeO3 thin film, dimensions: (|209)>,
          <EELSSpectrum, title: LaFeO3 thin film, dimensions: (|209)>,
          <EELSSpectrum, title: Strontium ferrite, dimensions: (|1024)>]
```



# Back to the demo...



### Load data downloaded from EELS Database

Core-loss and low-loss boron nitride spectra:

#### Set microscope parameters:

```
In [16]: s.set_microscope_parameters(beam_energy=100, convergence_angle=0.2, collection_angle=2.55)
```

#### Add Boron and Nitrogen to model:

```
In [17]: s.add_elements(('B', 'N'))
```






- Automatically adds PowerLaw component to fit background
- Automatically adds relevant edges in the energy range of the spectrum

```
In [18]: m = s.create_model(11=11)
         m.components
Out[18]:
                      Attribute Name
                                            Component Name
                                                                   Component Type
                                                                         PowerLaw
                            PowerLaw |
                                                  PowerLaw |
            1 |
                                 N_K |
                                                       N_K |
                                                                       EELSCLEdge
                                                       B_K |
                                                                       EELSCLEdge
                                 B_K |
```










# Perform least squares fitting to data

```
In [20]: m.smart_fit()
m.plot()
```






### Perform least-squares fitting with fine structure

- Requires local database of ionization cross-sections
- HyperSpy uses values calculated using Egerton's sigmak3 and sigmal3 routines
- Better ones come with DigitalMicrograph, and HyperSpy can use those files:
  - Not open-source, so cannot be distributed

```
In []: # To change to Gatan's GOS files:
    ### hs.preferences.EELS.eels_gos_files_path = <path to GOS tables>
In [21]: m.enable_fine_structure()
    m.smart_fit()
```











Next demo: <u>Processing TEM EDS data</u>

